

SMARTLINE
MEDICAL

Endoscopy Cabinet Training

iQ Series

Redefining Endoscopy Drying



ENDOSCOPY CABINET RANGE

iQ Series

Core Series



RotaScope Core
(not EN compliant)

Classic Series



SlidaScope
RotaScope



SlidaScope iQ
RotaScope iQ

OVERVIEW

Feature	Core	Classic	iQ Series
Integrated Air System Fully independent air system. IAS requires no infrastructure from the hospital and uses patented vacuum technology	✓	✓	✓
Plug and Play No complex install requirements	✓	✓	✓
Customisable Configuration (single door / Pass-Thru)	✓	✓	✓
EN16442 Compliant		✓	✓
Scope Data Monitoring & Intelligent Environmental Monitoring Touch Screen Display shows remaining scope storage time and environmental monitoring		✓	✓
Intelligent Interior LED Lighting		✓	✓
Remote Log In 24/7 service capability and reporting		✓	✓
Smart Storage (SlidaScope) Storage location based on preference		✓	✓
Channel Monitoring Individual scope lumen monitoring of humidity and airflow to prove scope dryness and safety			✓



* See Comparison Matrix

iQ Series

FEATURE BREAKDOWN

Channel Monitoring

World first technology

Validates drying of each scope channel in real time

Ensures patient safety

Ensures all channels are correctly connected

Smart Technology

Prioritises scope usage for user ergonomics

Intelligent lighting

Reminders for cleaning, servicing and PQ

Self monitoring environmental parameters

Intelligent Design



Fully Integrated Air System

No hospital infrastructure required


Small Footprint

Designed with looks and function in mind

CUSTOMIZATION



INTEGRATED AIR SYSTEM



- Fully independent air system
- Requires no infrastructure from the hospital
- Patented vacuum technology rules out possibility of cross contamination
- Air flow with individual channel drying and world first **individual channel monitoring**.

EN16442 Air Quality Certification

- Air Quality-moisture content:
 - Subclause 5.2.1.1.2 test 6.6.2 – Pass
- Air Quality-oil content:
 - Subclause 5.2.1.1.2 test 6.6.3 – Pass/optional
- Air Quality particulate contamination:
 - Subclause 5.2.2.4 test 6.6.1 & 6.11 – Pass
- Airborne microbial contamination:
 - Subclause 5.2.2.1 test Annex C – Pass



iQ Series


ACTIVATING TRACKING

Easy 3 step process at washer and cabinet

1. **WHO** – Operator barcode presented to station
2. **WHAT** – Scan Scope barcode
3. **WHERE** – Scope loaded or unloaded from device



AER
















CABINET



iQ Series

TRACKING AND TRACEABILITY

ENVIRONMENTAL	LOCATION	DRY STATUS	ENDOSCOPE NAME	REMAINING TIMES:	STORAGE	OoC
TEMP 	1		PCF180 2236524		167:44	2:43
VAC 	2		CF190 26665322		76:33	2:12
AIR 	3					
PUMP TEMP 	4		GIF190 2234572		122:42	2:48
F  R	5		PCF180 2236528		101:54	0:00
HUMIDITY 	6		GIF190 2557464		163:44	2:43
CLEAN 	7					
Clean Due 22/12/2020	8		CF190 29654344		165:45	0:15
CONTROL PANEL	9					

- Barcode tracking system
- Touch screen display shows remaining scope validity storage time in and out of cabinet
- Compatible with any AER model
- Paperless validation of scope drying

The Environmental Monitor measures and records the following cabinet conditions:

- Temperature
- Humidity
- Air (Cabinet Pressure)
- Vacuum
- Door Status
- IAS Temperature
- Routine Cleaning Status

iQ Series

TRACKING AND TRACEABILITY – CHANNEL MONITORING

WASHED TIME
12-Jul-22 - 14:49

WASHED BY
NA

STORED BY
TAYLOR

CHANNEL	DRY STATUS	SCOPE LOCATION: 3
C1		N/A
C2	⚠️	CHECK CONNECTION
C3		N/A
C4	☀️	DRY
C5	☁️🌧️	
C6	☀️	DRY
C7	☁️🌧️	

EXPIRE SCOPE /
FORCE REMOVAL

CH2 Connection

- Shows individual Channel Status
- Full scope details available
- Records errors and faults
- Paperless validation of scope drying

iQ Series

CONTROL SCREEN

CONTROL PANEL

HOME

HELP

SENSORS

ALARM
ACKNOWLEDGE

Power-Reset (Purple light),

AIR OFF

LOCK CABINET

REMOVE ALL SCOPES
FOR CLEANING

Version: 1.0



Data Capture

Theatre Use

- Patient



Bedside Clean

- Immediate flush and wipe down



Drying Cabinet



AER Washer

- Automated Endo Reprocessor
- Chemical Cleaning
- Sometimes Alcohol flush and air purge



Manual Decontamination

- Sink cleaning
- Leak testing



iQ Series

TICKET PRINTOUT


- Scopes can be presented to printer at any stage for instant status report
- Tickets can be customised based on site preference
- Printed onto stickers for easy physical record keeping and daybook procedures

SCOPE WASH RECORD

SCOPE TYPE: ENDOSCOPE
SCOPE NAME: scope 4
WASHER: AER 1
LOADED TIME: 11:46:36 Jan-06 20
LOADED BY: Richard2
REMOVAL TIME: 11:46:42 Jan-06 20
REMOVED BY: Richard2
MAX STORAGE: 168 HOURS
OOC EXPIRES: 14:46:52 Jan-06 20

SCOPE STORAGE RECORD

SCOPE TYPE: ENDOSCOPE
SCOPE NAME: scope 4
LOADED TO: ROTASCOPE 1:1
LOADED TIME: 11:46:58 Jan-06 20
LOADED BY: Door TAG
REMOVAL TIME: 11:47:09 Jan-06 20
REMOVED BY: Richard2
STORED FOR: 0.0 HOURS
OOC EXPIRES: 14:46:52 Jan-06 20



Site Estimation

Considerations are listed in the **F227** form. This must be completed for every cabinet job and given to Smartline before production commences to ensure job details are correct and assist in preparing the site for installation.

Site Visit and Consultation



```
graph TD; A[Site Visit and Consultation] --> B[Quote and Design - Preparation and Presentation]; B --> C[Customer Order]; C --> D[Confirmation of Requirements (F227) and Site preparation]; D --> E[Installation and Training];
```

Quote and Design - Preparation and Presentation

Customer Order

Confirmation of Requirements (F227) and Site preparation

Installation and Training

Data Equipment

G-Series Cabinet +

SYSTEM STARTER KIT

Limited functionality without adding more Washer Stations
Limited Qty Plaques included



LAPTOP KIT



PRINTER KIT



WASHER EXPANSION KIT



iQ and Classic Cabinet +

COMMUNICATION KIT



SMART STATION

Plus 2 Cables
One Station
does everything with 10m Radius



+ FREE OF CHARGE PER JOB
All Barcode labels



Communication Kit


- Wifi Router
- Laptop (optional)
- Additional Wifi Dongle
1x per site

Smart Station

- Printer
- Touch Screen
- Barcode Scanner
- Wall Mount
1+ per site (as required)

Barcodes

- Washer
- Operator
- Scopes



Case Study: Single Door Cabinets

Locations for:

- Cabinets
- Smart Stations
- Data Requirements
- Power Requirements

Considerations:

- Ceiling Height
- Services/Obstructions
- Access

EQUIPMENT TO BE QUOTED

Total Requirements


- 6x cabinets
- 1x Communication Kit
- 1x Smart Station
- Hose Kits?



EXISTING HEIGHT ADJUST
S/S SINK BENCH & OVER
BENCH SELF TO BE
RETAINED

EXISTING LAM CUPBDS &
OVER BENCH SHELVES
TO BE RETAINED

MIRIAN HOUSE, BLOCK M



Case Study: Pass Thru Cabinets

Locations for:

- Cabinets
- Smart Stations
- Data Requirements
- Power Requirements

Considerations:

- Ceiling Height
- Services/Obstructions
- Access
- Wall Cut Out

EQUIPMENT TO BE QUOTED

Total Requirements

- 4x cabinets (2x Single Door, 2x Pass thru)
- 1x Communication Kit
- 3x Smart Stations
- Hose Kits?




COMMUNICATION KIT



Case Study: Pass Thru Cabinets

- Data Requirements
- Power Requirements
- Wall Cut Out



Hose Kits (Scope Connectors)

Hose kits are included FOC with each cabinet purchase for standard scopes.

- Endoscopes from different brands and of different types will have ports that are unique and require specific connectors

Some specialty scopes will require quoting (Not FOC)

- EUS/EBUS/Duodenoscope/etc

Please ensure you get a scope fleet list from your customer with scope models and quantities so this can be assessed early.



SMARTLINE
MEDICAL

iQ Series

Redefining Endoscopy Drying

