

F227 - Support Document

Reference File for completing an F227 Form



INTRODUCTION

This document is intended to assist with gathering required information and details when accepting a purchase of a EN16442 Compliant endoscope cabinet.

This may also be used to assist in accurate quoting at the beginning of the customer engagement and can be shared with the customer when required.

The F227 Document is to be completed by distributors and is intended to assist in cleanhandover of the EN16442 compliant cabinet range.

SITE DETAILS

Hospital: Site the cabinet is intended for

Purchase Order: Reference number

ROOM INFORMATION

Ceiling Height to be considered: Rotascope needs a clearance height of 2500mm or more.

Revision of wall opening drawings need to be supplied to the customer for pass through cabinets. These can be found in Distributer Portal.

Pass Through wall thickness: If the wall thickness is more than 300mm a custom manufactured solution is possible and needs to be arranged prior to install. A custom version will have a longer lead-time and will be at additional cost.

CABINET INFORMATION

Quick reference to the size of the job, with more detail about each cabinet to follow in cabinet details part of the F227.

Quantity of cabinets to be recorded.

Quantity of stickers to be recorded for sides of cabinet.

Photo of barcode type if there are existing ones already in use.

SCOPE INFORMATION

Capture brand, model and quantities of each scope intended for use in cabinet.

Smartline will determine from the information provided which kitsets are needed.

DATA COMPONENT

Smart station quantities and communication kits required.

Pass Through cabinets require 2 x smart stations each and communication kits only required for IQ & Classic Cabinets.

Any additional information needed for data set-up.

CABINET DETAILS

Each cabinet is to be itemized to accommodate for cabinet models (the cabinets codes are listed on the F227 form).

Specify sticker option if cabinet is an IQ Series Cabinet. Classic & Core Cabinets have standardized stickers, but the IQ Series has 4 different options to choose from (see Appendix A of F227 form).

Door hinge side to be added to this area, if cabinet is a pass through both front and rear hinge sides need to added.

Include any additional information relating to each cabinet here.