

SlidaScope
(09.DC-B-SS)
2018 - 2021 MODEL



SlidaScope Classic
(SLIDA-S)
2021- onwards MODEL



SlidaScope iQ
(SLIDA-IQ-S)
2021- onwards MODEL



Size	600mmW x 600mmL x 1820mmH (2140mmH with SAB)	600mmW x 600mmL x 1880mmH	600mmW x 600mmL x 1880mmH
Air System	Shared Air Box (SAB)	Inbuilt Integrated Air System (IAS)	Inbuilt Integrated Air System (IAS)
Servicing	Alternating minor and major service	One annual service: 2x pumps, 1x HEPA filter, 1x intake filter	One annual service: 2x pumps, 1x HEPA filter, 1x intake filter
Pass Through	Flashing kit required for rear door	Cabinet mounted rear door	Cabinet mounted rear door
Channel Monitoring	No	No	Yes
Data System	RFID	Barcode	Barcode
Hose Kits Included	Olympus (standard)	All available (Olympus standard scope kits included FOC)	All available (Olympus standard scope kits included FOC)
Price	Refer Price List (CABINET + SAB)	Refer Price List	Refer Price List