

SETTING NEW STANDARDS



redefining endoscopy drying



INTRODUCING

CHANNEL MONITORING

World first technology

Validates drying of each scope channel in real time

Confirms all channels are correctly connected

Ensures Patient Safety



CABINET FEATURES

- **Integrated Air System**, no hospital infrastructure required
- **Smart Technology including unique Channel Monitoring:** validating scope drying
- **Customisable** and available in Single Door, Pass Through or X-Ray Pass Through
- **Smart environmental monitoring**
- **EN16442 compliant**, approved for 31 day storage (7 days in AUS/NZ)
- **Smart LED interior lighting**, colour coded
- **Automated cleaning**, service, and PQ reminders
- **Touch screen display** shows remaining scope validity storage time using barcode
- **Compatible** with any AER model

ROTASCOPE

- **Intelligent Endoscope Storage Cabinet:** Stores up to 9 flexible endoscopes vertically
- **Small footprint.** Dimensions: 600mmW x 600mmL x 2363mmH

SLIDASCOPE

- **Intelligent Endoscope Storage Cabinet:** Stores up to 8 flexible endoscopes horizontally in trays
- **Small footprint.** Dimensions: 600mmW x 600mmL x 1862mmH



watch the video here



GOOD DESIGN AWARD WINNER



+61 7 5478 9977 • www.smartlinemedical.com • 55 Cordwell Road, Yandina QLD 4561

SMARTLINE MEDICAL